



Flash Matrix

All Results in One System



Random
Access



Adjustable
Monitor



Mobile & Small
Footprint



Flash Matrix

Model: M-16 / M-32



Product Features

High Throughput: 16~64 Samples

Flexible: Random Access

Fast: 1 min Sample Prep, 30 mins Sample to Result

Simple: All Instruments Connects with One Monitor, Less Operation & Maintenance

Easy: All Results in One System, Operator Friendly

Technical Specifications

Product Name	Flash Matrix	
Model	M-16	M-32
Instruments Equipped	4 in 1	8 in 1
Test Throughput	16-32 (Max.64)	
Instruments Model	Flash20 (Customization is available)	
Fluorescence Channel	FAM/SYBR Green(470nm), HEX/VIC (525nm), ROX/TEXAS RED (580nm), CY5(635nm)	
Detection Probe Temperature	High Sensitivity Photodetector	
Control Mode Data	Independent Temperature Control	
Communication Interface	USB	
Power Requirements	220V	
Size	L74*W60*H116 cm	L74*W60*H147 cm
Weight	30 kg	60 kg

Coyote Bioscience Co., Ltd.

H.O.: 5/F, Chuangyezhong Road #36, Haidian District,
Beijing 100085, P. R. China.

Tel: +86 10 8071 8308 ext.8070

Email: sales@coyotebio.com



www.coyotebio.com